paediSHUNTASSISTANT® / SHUNTASSISTANT®

Ⓓ Patientenhandbuch | ⒸⒺ ⒺⒶ Patient Manual | ⒻⒼ ⒼⒶ Manuel du patient
ⒺⒺ ⒺⒶ Manual para el paciente | ⒺⒺ ⒺⒶ Manuale per il paziente

CAUTION

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THE COMPANY

Christoph Miethke GmbH & Co. KG is a Brandenburg-based company that develops, manufactures and markets innovative neurosurgical implants for the treatment of hydrocephalus. In this, we work in successful partnerships with numerous hospitals worldwide.

The purpose of this booklet is to provide you and your family with some understanding of the treatment of hydrocephalus. The successful treatment of this condition has only been possible since the 1950s. In a dramatic race against time to save the life of his son, Casey, who suffered from hydrocephalus, a technician named John D. Holter developed, in only a few weeks, a novel silicone valve. Despite the fact that, since its first implantation in March 1956, this valve has proven to be clinically effective and a giant step in the treatment of this condition, there are many patients today who experience considerable problems with hydrocephalus valve systems.

Christoph Miethke picked up the knowledge gained in 50 years of valve treatment and developed a new generation of highprecision valves made of the metal titanium. For the first time, there are valve systems available that consistently take into account the physical conditions of brain fluid drainage and can thus maintain a physiological brain pressure, independent of the body position of the patient.

Fig. 1: Anatomic sketch of the cranium (inner cover page)

- 1 skull pan
- 2 brain
- 3 cerebrospinal fluid
- 4 lateral ventricle
- 5 third ventricle
- 6 fourth ventricle

BASIC ANATOMY

The human brain (Fig. 1) is surrounded by a special fluid known as cerebrospinal fluid (CSF). Cerebrospinal fluid is produced in several chambers, so-called ventricles, that are found within the brain. The channels, by which the ventricles are interconnected, constitute a complex drainage system. The fluid in the brain circulates through these ventricles and eventually flows into the venous blood. The function of this fluid is to protect the brain from mechanical damage. The CSF also regulates the internal brain pressure (intracranial pressure, ICP), keeps the brain tissue moist and transports the products of metabolism.

CLINICAL PICTURE OF THE CONDITION

In healthy humans, a balance exists between the production and resorption of cerebrospinal fluid. In infants, approx. 100 ml of this fluid is produced every day; in small children, the daily production is approx. 250 ml, in adults approx. 500 ml. If the amount of fluid produced exceeds the amount resorbed, the ventricles expand, leading to the condition known as hydrocephalus (Fig. 2). The term hydrocephalus refers to the continuous increase of the volume of "water" (hydro) in the "head" (cephalus). This condition is often observed at birth (congenital hydrocephalus), but it can also develop later in life, e.g., as the result of inflammation, hemorrhage or severe head injury, or after brain surgery. Such cases are referred to as acquired hydrocephalus.

A further distinction is made between obstructive hydrocephalus and communicating hydrocephalus. In obstructive hydrocephalus, the links between the ventricles of the brain are interrupted so that the ventricles cannot "communicate" with each other. Cases in which the ventricles are interlinked through open channels, but resorption of cerebrospinal fluid is impaired, are diagnosed as communicating hydrocephalus.


Fig. 2: Ventricle size
a) normal, b) hydrocephalus

CLINICAL SYMPTOMS OF THE CONDITION

In infants, the cranial bones have not adhered solidly yet. The increasing volume of cerebrospinal fluid causes the head to increase in circumference while, at the same time, brain tissue disintegrates. From the age of about 2, the hardened skull prevents any growth of the head's circumference. In that situation, the increase in fluid volume leads to a massive pressure rise, resulting in the expansion of the brain ventricles and the compression of the brain itself. The consequence for infants and adults

can be irreversible brain damage. Symptoms (depending on the severity of the disorder) include nausea, headache, vomiting, impaired coordination, drowsiness and, in the end, unconsciousness.

DIAGNOSIS OF THE CONDITION

Doctors have a variety of ways at their disposal to diagnose hydrocephalus. The ventricle size is measured through imaging procedures (e.g. computerized tomography, ultrasound or NMR-tomography).

Computerized tomography (CT)

This quick and painless diagnostic procedure produces X-ray images of different layers of the head.

Nuclear Magnetic Resonance (NMR) tomography

This painless electromagnetic imaging process produces images of very fine layers of the head. It is also known as NMR, MRT, or MRI scanning.

Ultrasound

This procedure, in which the interior of the head is examined through the open fontanel, can only be applied to small children.

Another way of diagnosing hydrocephalus is through pressure measurements showing an increased brain pressure. The circulation of cerebrospinal fluid is investigated through examinations with contrast agents.

METHODS OF TREATMENT

For all the efforts to find therapeutic alternatives to valve implantation (e. g. through pharmaceutical treatment or, most recently, by minimally invasive surgery), there is currently no alternative, in most cases, to the implantation of a drainage system, referred to as a shunt. The operation is usually neither risky nor difficult. The drainage system (Fig. 3) is comprised of catheters that drain off the cerebrospinal fluid and a valve that regulates intracranial pressure. In many cases, a reservoir is implanted between the ventricle catheter and the valve. The doctor uses this reservoir for checking the patency of the valve, to siphon off CSF-samples or to inject medicines.

A distinction is made between two types of drainage: ventriculo-peritoneal (from the head to the abdominal cavity) and ventriculo-atrial (from the

head to the right atrium of the heart). In special cases, there is also the option of implanting a lumboperitoneal shunt (from the spine canal into the abdominal cavity).

- 1 right atrium
- 2 heart catheter (atrial catheter)
- 3 valve
- 4 reservoir
- 5 ventricular catheter
- 6 ventricles
- 7 abdominal catheter (peritoneal catheter)
- 8 abdominal cavity


Fig. 3: Drainage systems for hydrocephalus patients
a) ventriculo-atrial, b) ventriculo-peritoneal

THERAPY COMPLICATIONS

A number of different complications can arise from shunt therapy for hydrocephalus, including infections, obstruction of the drainage system and undesirably excessive drainage of cerebrospinal fluid. This can cause the valve system to induce obstruction or increased drainage. The *paediSHUNTASSISTANT/ SHUNTASSISTANT* in combination with a simple or adjustable valve, avoids complications that might occur as the result of increased drainage of cerebrospinal fluid. This is possible because, in addition to the valve, the *paediSHUNTASSISTANT/ SHUNTASSISTANT* also regulates intracranial pressure in vertical and sitting position. This allows for the achievement of a state that closely approximates the pressure ratios in normal individuals. To give you an understanding why your physician decided for the *paediSHUNTASSISTANT/ SHUNTASSISTANT*, see chapter Physics Background.

AFTER THE OPERATION

As a rule, the everyday activities of patients with shunt implants are not restricted. However, patients should abstain from major physical exertion (e. g. hard physical work, strenuous sports). Hydrocephalus patients who experience headache, dizziness, unnatural gait or similar symptoms should consult a physician without delay. Apart from that, we recommend medical check-ups at regular intervals. The patient should avoid knocks or pressure on the valve and catheters. The valve has been designed to be resistant against magnetic fields.

PHYSICS BACKGROUND

In the following chapter we describe the pressure conditions relevant for hydrocephalus drainage. The ventricle pressure and the pressure in the abdominal cavity are represented by water levels. In a healthy human, the ventricular pressure (water level in the ventricle container) is positive (slightly above 0) in the horizontal position and negative (slightly below 0) in the upright (vertical) position (Fig. 4).

- 1 ventricle container
- 2 ventricles


Fig. 4a: Ventricle pressure in a healthy human in horizontal position


Fig. 4b: Ventricle pressure in a healthy human in vertical position

- 1 ventricle container
- 2 expanded ventricles


Fig. 5a: Ventricle pressure in a hydrocephalus patient in horizontal position


Fig. 5b: Ventricle pressure in a hydrocephalus patient in vertical position

Now there is an urgent need to lower the intracranial pressure and keep it within normal limits, independent of the body position. Finally, excessive cerebrospinal fluid is drained into the abdominal cavity. Due to changes in the patient's body position, the drainage system is subject to incessant, considerable physical changes. Fig. 6 shows the effects on the intracranial pressure when a tube is implanted, although there has been no valve integrated in the drainage system yet.

For simplicity, the abdominal pressure as well as the ventricles can be regarded as open vessels, which are now connected by a tube. As long as the patient is lying down (head and abdomen at the same height) and no valve is integrated in the drainage system, both water levels are at the same height, too: It is a system of communicating vessels. In a simplified picture, the abdomen can be regarded as an overflow vessel. Even if more fluid is filled into the ventricle container, the water level in it will remain at the same height, because the fluid is instantly drained into the abdomen. When the patient stands up, the ventricles are at a significantly higher level than the abdomen. In this case, the

fluid is drained through the tube until both water levels are at the same height. This means, however, that the ventricle container is emptied completely. Since the ventricles do not have rigid walls, this emptying leads to a contraction of the ventricles. This, in turn, can result in the above mentioned blockage of the drainage system. The cerebrospinal fluid is sucked out of the ventricles and the brain suffers deformation. And when the brain contracts, the resulting open space between brain and skull bone can fill up with water or blood (Fig. 6).

- 1 ventricle container
- 2 ventricle
- 3 drainage tube
- 4 abdominal cavity
- 5 abdominal container
- 6 accumulated water or blood
- 7 contracted ventricles



Fig. 6: Ventricle drainage without a valve
a) horizontal, b) vertical

A conventional valve integrated into the drainage system causes a rise of the water level in the ventricle container, by exactly the opening pressure of the valve. Now, the two containers "communicate" only when the valve is open. When the patient stands up, cerebrospinal fluid is drained off until the height difference between the two containers for the horizontal position is reached. However, the opening pressure of the valve, which was adjusted for the horizontal position, is considerably lower than the pressure corresponding to the height difference between the ventricles and

the abdomen. Hence, the ventricles will still be drained empty, resulting in the above mentioned problems (Fig. 7).

- x opening pressure in horizontal position
- 1 ventricle container
- 2 ventricle
- 3 drainage tube
- 4 abdominal cavity
- 5 abdominal container
- 6 accumulated water or blood
- 7 contracted ventricles
- 8 valve



Fig. 7: Ventricle drainage with conventional valve
a) horizontal, b) vertical

The simplified diagram clearly demonstrates how crucial it is that a valve is implanted whose opening pressure is considerably higher for standing (and sitting) position (depending on the distance between the brain and abdomen) than for horizontal position. If a simple or adjustable valve is nonetheless implanted that only functions adequately in horizontal position, the shunt system must be equipped with a supplementary valve for standing (and sitting) positions.

The *paediSHUNTASSISTANT*/ *SHUNTASSISTANT* is such a supplementary valve. It functions exclusively when the patient is standing (or sitting) and, in combination with the differential pressure valve, ensures that the intracranial pressure will be exactly what the patient needs in every body position. The problems and complications caused by undesirably high levels of drainage are thus avoided (Fig. 8).

- x opening pressure in horizontal position
 1 ventricle container
 2 ventricles
 3 drainage tube
 4 abdominal cavity
 5 abdominal container
 6 paediSHUNTASSISTANT/ SHUNTASSISTANT


Fig. 8: Ventricle drainage with paediSHUNTASSISTANT/ SHUNTASSISTANT
 a) horizontal, b) vertical

TECHNICAL SPECIFICATIONS OF THE VALVE

The valve is made of high quality materials that have been tested and standardized for use in implants. Materials used: titanium, tantalum, sapphire

Fig. 9: Schematic cross section of the paediSHUNTASSISTANT/ SHUNTASSISTANT (inner cover page)

- 1 titanium casing
 2 tantalum ball
 3 sapphire ball
 4 ball seat
 5 inlet port
 6 outlet port

The paediSHUNTASSISTANT/ SHUNTASSISTANT is implanted in conjunction with a simple differential pressure valve. It functions only in standing and sitting positions, but not in horizontal position (the differential pressure valve takes over in this position). The weight force of the tantalum ball compensates for the hydrostatic pressure difference that engenders undesirably high levels of cerebrospinal fluid drainage when the patient stands up, if a simple differential pressure valve is used. The paediSHUNTASSISTANT/ SHUNTASSISTANT opens if

critical intracranial pressure is exceeded while the patient is in standing position. A stable housing helps to reduce the effects of both external and internal factors (subcutaneous pressure, fluid) on valve function to a bare minimum.


Fig. 10 paediSHUNTASSISTANT/ SHUNTASSISTANT, Scale 1:1

Ideally, in adults the SHUNTASSISTANT is implanted in the thoracic cavity in accordance with anatomical relationships and in infants behind the ear in view of the problem of growth. The valve can be easily and safely implanted beneath the skin without limiting the scope of activities or movements performed hitherto by the patient.

Warning note: The shunt system may comprise a reservoir that can be pumped. Since frequent pumping may lead to overdrainage and thus to very unfavourable pressure conditions, such pumping should only be carried out by the physician.

PATIENT ID

Each paediSHUNTASSISTANT/ SHUNTASSISTANT is supplied with an individual patient ID. The ID is filled out by the attending physician and contains important information for the follow-up examinations.

A BRIEF PATIENT GLOSSARY

Anatomy

A guide to the structure of body components

Arachnoid

Connective tissue in the brain that lies between the dura mater and the pia mater

Catheter

Tube

Cerebrospinal fluid

Watery spinal fluid in the brain

Communicating vessels

Vessels that are connected by a channel

Computed tomography (CT)

Imaging technique whereby „slices“ of the body are recorded with an X-ray scanner

Drainage

Drainage of accumulated fluid

Dura mater

The hardest component of the meninx

Fluid component

Cerebrospinal fluid component

Fontanel

Connective tissue opening in a young infant's skull that later ossifies

Implant

Substance that is placed in the human body to replace a particular function for a limited period of time or for the rest of the patient's life

Lumbar puncture

Puncture of the spinal channel at the lower spine

Lumboperitoneal drainage

Drainage of cerebrospinal fluid from the ventricle of the brain, by way of the region of the lumbar vertebrae in the abdominal cavity

Meninges

Membrane found in the brain and spine

Meningitis

Inflammation of the meninx

Minimally invasive

Minimally infiltrating

Overdrainage

Undesirable outward flow of cerebrospinal fluid

Peritoneum

Membrane that covers the pelvic and abdominal cavities

Pia mater

Component of the soft meninx containing blood vessels

Piaarachnoid

Soft component of the meninx that is divided into the arachnoidea and pia mater

Puncture

Insertion of a hollow needle or a trocar into a vessel for the purpose of removing fluid

Resorption

Suctioning or removal of material through skin, mucosa or tissue

Shunt

A passage between two channels – here a catheter drainage system with an integrated valve

Spinal column

Element of the central nervous system located within the vertebral channel

Subcutaneous pressure

Pressure beneath the skin

Subdural hematoma

An accumulation of blood between the brain and cranium

Ventricle of the brain

Intracranial space containing cerebrospinal fluid

Ventricular peritoneal drainage

Drainage of cerebrospinal fluid from the ventricle of the brain directly into the abdominal cavity (abdominal catheter)

FOLLOW-UP EXAMINATIONS

A follow-up examination must be carried out in all cases.

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany