*proSA*<sup>®</sup>

Ⓓ Patientenhandbuch | ⒼⒹ ⒺⒶ Patient Manual | ⒻⒼ Manuel du patient  
ⒺⒼ Manual para el paciente | ⒾⒶ Manuale per il paziente

CAUTION

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## THE COMPANY

Christoph Miethke GmbH & Co. KG is a company based in Berlin-Brandenburg region that develops, manufactures and markets innovative neurosurgical implants for the treatment of hydrocephalus. In the course of our work we established successful partnerships with numerous hospitals worldwide.

The purpose of this booklet is to provide you and your family with some understanding of the treatment of hydrocephalus. The successful treatment of this condition has only been possible since the 1950s. In a dramatic race against time to save the life of his son, Casey, who suffered from hydrocephalus, a technician named John D. Holter developed, in only a few weeks, a novel silicone valve. Despite the fact that, since its first implantation in March 1956, this valve has proven to be clinically effective and a giant step in the treatment of this condition, there are many patients today who experience considerable problems with hydrocephalus valve systems.

Christoph Miethke utilised the knowledge gained in 50 years of valve treatment and developed a new generation of highprecision valves made of the metal titanium. For the first time, there are valve systems available that consistently take into account the physical conditions of brain fluid drainage and can thus maintain a physiological brain pressure, independent of the body position of the patient.

Fig. 1: Anatomic sketch of the cranium  
(see inner cover page)

- 1) skull
- 2) brain
- 3) cerebrospinal fluid
- 4) lateral ventricle
- 5) third ventricle
- 6) fourth ventricle

## BASIC ANATOMY

The human brain (fig. 1) is surrounded by a special substance known as cerebrospinal fluid (CSF). Cerebrospinal fluid is produced in several chambers, so-called ventricles, that are found within the brain. The channels, by which the ventricles are interconnected, constitute a complex drainage system. The fluid in the brain circulates through these ventricles and eventually flows into the venous blood. The function of this fluid is to protect the brain from mechanical damage. The CSF also regulates the internal brain pressure (intracranial

pressure, ICP), keeps the brain tissue moist and transports the products of metabolism.

## CLINICAL PICTURE OF THE CONDITION

In healthy people, a balance exists between the production and resorption of cerebrospinal fluid. In infants, approx. 100 ml of this fluid is produced every day; in small children, the daily production is approx. 250 ml, in grown-ups approx. 500 ml. If the amount of fluid produced exceeds the amount resorbed, the ventricles expand, leading to the condition known as hydrocephalus (fig. 2). The term hydrocephalus refers to the continuous increase of the volume of "water" (hydro) in the "head" (cephalus). This condition is often observed at birth (congenital hydrocephalus), but it can also develop later in life, e.g., as the result of inflammation, hemorrhage or severe head injury, or after brain surgery. Such cases are referred to as acquired hydrocephalus.


Fig. 2: Ventricle size  
a) normal, b) hydrocephalus

A further distinction is made between obstructive hydrocephalus and communicating hydrocephalus. In obstructive hydrocephalus, the links between the ventricles of the brain are interrupted so that the ventricles cannot "communicate" with each other. Cases in which the ventricles are interlinked through open channels, but resorption of cerebrospinal fluid is impaired, are diagnosed as communicating hydrocephalus.

## CLINICAL SYMPTOMS OF THE CONDITION

In infants, the cranial bones have not yet grown together solidly. The increasing volume of cerebrospinal fluid causes the head to increase in circumference while, at the same time, brain tissue disin-

tegrates. From the age of about 2, the hardened skull prevents any growth of the head's circumference. In this case, the increase in fluid volume leads to a massive pressure increase, resulting in the expansion of the brain ventricles and the compression of the brain itself. The consequence for infants and adults can be irreversible brain damage. Symptoms (depending on the severity of the disorder) include nausea, headache, vomiting, impaired coordination, drowsiness and, in the end, unconsciousness.

## DIAGNOSIS OF THE CONDITION

Doctors have a variety of methods at their disposal to diagnose hydrocephalus. The ventricle size is measured through imaging procedures (e.g. computerized tomography, ultrasound or NMR-tomography).

### Computerized tomography (CT)

This quick and painless diagnostic procedure produces X-ray images of different layers of the head.

### Nuclear Magnetic Resonance (NMR) tomography

This painless electromagnetic imaging process produces images of very fine layers of the head. It is also known as NMR, MRT, or MRI scanning.

### Ultrasound

This procedure, in which the interior of the head is examined through the open fontanel, can only be applied to small children.

Another way of diagnosing hydrocephalus is through pressure measurements showing an increased brain pressure. The circulation of cerebrospinal fluid is investigated through examinations with contrast agents.

## METHODS OF TREATMENT

For all the efforts to find therapeutic alternatives to valve implantation (e. g. through pharmaceutical treatment or, most recently, by minimally invasive surgery), there is currently no alternative, in most cases, to the implantation of a drainage system, referred to as a shunt.

- 1 right atrium
- 2 heart catheter (atrial catheter)
- 3 valve
- 4 reservoir
- 5 ventricular catheter
- 6 ventricles
- 7 abdominal catheter (peritoneal catheter)


Fig. 3: Drainage systems for hydrocephalus patients  
a) ventriculo-atrial, b) ventriculo-peritoneal

## THERAPY COMPLICATIONS

In certain instances, complications can result from the treatment of hydrocephalus with a shunt system. As with any surgical procedure, an infection can develop. Unfortunately, problems can also arise that are either directly or indirectly linked with the implanted valve system. These complications include blockage of the drainage system, an unavoidable adaptation of the system to the growth of the child, and an unintended increase in drainage of cerebrospinal fluid. In order to help you understand why your physician decided to use *proSA*, the basic physical principles underlying a typical drainage system are explained in the section entitled „Physics background“.

## AFTER THE OPERATION

As a rule, the everyday activities of patients with shunt implants are not restricted. However, patients should abstain from major physical exertion (e. g. hard physical work, strenuous sports). Hydrocephalus patients who experience headache, dizziness, unnatural gait or similar symptoms should consult a physician without delay. Apart from that, we recommend medical check-ups at regular intervals.

The patient should avoid knocks or pressure on the valve and catheters. The valve has been designed to be resistant against magnetic fields.

## PHYSICS BACKGROUND

The following chapter describes the pressure conditions relevant for hydrocephalus drainage. The ventricle pressure and the pressure in the abdominal cavity are represented by water levels. In a healthy human, the ventricular pressure (water level in the ventricle container) is positive (slightly above 0) in the horizontal position and negative (slightly below 0) in the upright (vertical) position (see fig. 4)

1 ventricle container  
2 ventricles


Abb. 4a: Ventricle pressure in a healthy human in horizontal position


Fig. 4b: Ventricle pressure in a healthy human in vertical position

In hydrocephalus patients, the ventricular pressure is always increased (water level in the ventricle container far above 0) regardless of the body position. The ventricles are expanded (see fig. 5).

1 ventricle container  
2 ventricles


Fig. 5a: Ventricle pressure in a hydrocephalus patient in horizontal position


Fig. 5b: Ventricle pressure in a hydrocephalus patient in vertical position

Now there is an urgent need to lower the intracranial pressure und keep it within normal limits, regardless of the body position. Finally, excessive cerebrospinal fluid is drained into the abdominal cavity.

Fig. 6 shows the effects on the intracranial pressure when a tube is implanted, although there has been no valve integrated in the drainage system yet. For simplicity, the abdominal pressure as well as the ventricles can be regarded as open vessels, which are now connected by a tube. As long as the patient is lying down (head and abdomen at the same height) and no valve is integrated in the drainage system, both water levels are at the same height, too: It is a system of communicating vessels. In a simplified picture, the abdomen can be regarded as an overflow vessel. Even if more fluid is filled into the ventricle container, the water level in it will remain at the same height, because the fluid is instantly drained into the abdomen.

When the patient stands up, the ventricles are at a significantly higher level than the abdomen. In this case, the fluid is drained through the tube until both water levels are at the same height. This means, however, that the ventricle container is emptied

completely. Since the ventricles do not have rigid walls, this drainage leads to a contraction of the ventricles.

- 1 ventricle container
- 2 ventricle
- 3 drainage tube
- 4 abdominal cavity
- 5 abdominal container


Fig. 6: Ventricle drainage without a valve  
a) horizontal, b) vertical

A conventional valve integrated into the drainage system causes a rise of the water level in the ventricle container, by exactly the opening pressure of the valve. Now, the two containers “communicate” only when the valve is open. When the patient stands up, cerebrospinal fluid is drained off until the height difference between the two containers for the horizontal position is reached. However, the opening pressure of the valve, which was adjusted for the horizontal position, is considerably lower than the pressure corresponding to the height difference between the ventricles and the abdomen. Hence, the ventricles will still be drained empty, resulting in the above mentioned problems (Fig. 7).


Fig. 7: Ventricle drainage with conventional valve  
a) horizontal, b) vertical,  $x$  = opening pressure

The simplified diagram clearly demonstrates how crucial it is that a valve is implanted whose opening pressure is considerably higher for standing (and sitting) position (depending on the distance between the brain and stomach) than for horizontal position. If a simple or adjustable valve is nonetheless implanted that only functions adequately in horizontal position, the shunt system must be equipped with a supplementary valve for standing (and sitting) positions.

The *proSA* is such a supplementary valve. It functions exclusively when the patient is standing (or sitting) and, in combination with the differential pressure valve, ensures that the intracranial pressure will be exactly what the patient needs in every body position. The problems and complications caused by undesirably high levels of drainage are thus avoided (Fig. 8).


Fig. 8: Ventricle drainage with proSA valve  
a) horizontal, b) vertical

## TECHNICAL SPECIFICATIONS OF THE VALVE

The *proSA* is a combination of an adjustable gravitational unit and a differential pressure unit. In certain cases, the adjustable unit can be implanted on its own. The dimensions are shown in Fig. 9.


Fig. 9: *proSA*, adjustable gravitational unit with  
(a) differential pressure unit or with  
(b) adjustable differential pressure unit

The closer the upper body of the patient is to the vertical position, the higher is the total opening pressure of the *proSA*.

This opening pressure can be changed through the skin even after the valve has been implanted, without necessitating any further surgical intervention. The *proSA* is made only of high-grade materials that have been tested and standardized for use in implants. The main material is titanium. The robust housing of the valve ensures that any influences (e. g. external pressure) on the functioning of the valve are reduced to an insignificant minimum. This guarantees a high level of functional safety and a long service life.

Fig. 10: Schematic cross section of the *proSA*  
(see inner coverpage)

- 1 titanium casing
- 2 rotor
- 3 bow spring
- 4 weight
- 5 sapphire ball

**Warning note:** The shunt system may comprise a reservoir that can be pumped. Frequent pumping can lead to overdrainage and thus pressure conditions outside of the normal physiological range. As such pumping should only be carried out by the physician.

## PATIENT ID

Each *proSA* is supplied with an individual patient ID. The ID is filled out by the attending physician and contains important information for the follow-up examinations.

## A BRIEF PATIENT GLOSSARY

### Anatomy

A guide to the structure of body components

### Arachnoid

Connective tissue in the brain that lies between the dura mater and the pia mater

### Catheter

Tube

### Cerebrospinal fluid

Watery spinal fluid in the brain

### Communicating vessels

Vessels that are connected by a channel

### Computed tomography (CT)

Imaging technique whereby „slices“ of the body are recorded with an X-ray scanner

### Drainage

Drainage of accumulated fluid

### Dura mater

The hardest component of the meninx

### Fluid component

Cerebrospinal fluid component

### Fontanel

Connective tissue opening in a young infant's skull that later ossifies

Implant Substance that is placed in the human body to replace a particular function for a limited period of time or for the rest of the patient's life

### Lumbar puncture

Puncture of the spinal channel at the lower spine

### Lumboperitoneal drainage

Drainage of cerebrospinal fluid from the ventricle of the brain, by way of the region of the lumbar vertebrae in the abdominal cavity

### Meninges

Membrane found in the brain and spine

### Meningitis

nflammation of the meninx

### Minimally invasive

Minimally infiltrating

### Overdrainage

Undesirable outward flow of cerebrospinal fluid

### Peritoneum

Membrane that covers the pelvic and abdominal cavities

### Pia mater

Component of the soft meninx containing blood vessels

### Piaarachnoid

Soft component of the meninx that is divided into the arachnoidea and pia mater

### Puncture

Insertion of a hollow needle or a trocar into a vessel for the purpose of removing fluid

### Resorption

Suctioning or removal of material through skin, mucosa, or tissue

### Shunt

A passage between two channels – here a catheter drainage system with an integrated valve

### Spinal column

Element of the central nerous system located within the vertebral channel

### Subcutaneous pressure

Pressure beneath the skin

### Subdural hematoma

An accumulation of blood between the brain and cranium

### Ventricle of the brain

Intracranial space containing cerebrospinal fluid

### Ventricular peritoneal drainage

Drainage of cerebrospinal fluid from the ventricle of the brain directly into the abdominal cavity (abdominal catheter)

FOLLOW-UP EXAMINATIONS

A follow-up examination must be carried out in all cases.

Date	Treatment

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Christoph Miethke GmbH & Co. KG | Ulanenweg 2 | 14469 Potsdam | Germany
Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany